

Plasma Metagenomic Sequencing & BALF Metagenomic Sequencing

**Leigh Furman, PharmD - Sr. Medical Science Liaison, Karius
APIC - Puget Sound Fall Conference - September 18th, 2025**

Background

Advantages and Limitations of Various Diagnostic Tools

	Culture	Antigen Tests	Serology	PCR
Advantages	- Traditional “gold standard” - Can isolate colonies - Can perform susceptibility	- Rapid - Inexpensive	May detect non-culturable organisms	- Relatively higher sensitivity - Relatively higher specificity - May have shorter TAT (in-house v. send-out) - Multiplex PCR panels 4 to > 25 pathogens - Some are FDA-cleared
Limitations	- Can have prolonged TAT - Limited sensitivity - Limited species resolution - May require invasive sample	- Variable sensitivity - Variable specificity (subject to false-positives, or cross-reactivity between pathogens)	- Lower sensitivity in IC hosts - Subject to pathogen cross-reactivity - Weeks to months for antibody production (false negative in early infection) - Prone to false-positives and/or persistent positivity	- Fewer pathogens detected - Require a-priori hypothesis for suspected organism - May require send-out

TAT = turnaround time; PCR = polymerase chain reaction; IC = immunocompromised

Plasma Metagenomic Sequencing (mNGS) Test Overview

Liquid biopsy for infectious diseases

Plasma mNGS is a blood-based test using metagenomic next-generation sequencing (mNGS) of **microbial cell-free DNA (mcfDNA)**.



1 blood sample
(specifically,
plasma)




1 day
processing




1000+ detectable
pathogens

Microbial cell-free DNA (mcfDNA)




1. Blauwkamp, T. A. et al. *Nat. Microbiol.* 4, 663–674 (2019)
2. Hu, P., Liang, D., Chen, Y. et al. *J. TranslMed.* 17, 124 (2019)
3. Volik, S. et al. *Mol. Cancer Res.* 14, 898–908 (2016)

Plasma mNGS Provides Minimally-invasive Detection of Pathogens Associated with Deep-seated and Difficult-to-diagnose Systemic Infections



Plasma Metagenomic Sequencing – Test Process



The Plasma Metagenomic Sequencing Test Can Detect Common AMR Markers for 4 Classes of Antimicrobials Across 18 Bacterial Pathogens

Gram-Negative Bacteria	AMR Marker						
	SCCmec	mecA	mecC	vanA	vanB	CTX-M	KPC
Carbapenem resistant and extended spectrum beta-lactamase (ESBL) producing bacteria							
<i>E. cloacae</i>						+	+
<i>E. coli</i>						+	+
<i>K. aerogenes</i>						+	+
<i>K. pneumoniae</i>						+	+
<i>K. oxytoca</i>						+	+
<i>P. mirabilis</i>						+	+
<i>P. vulgaris</i>						+	+
<i>S. bongori</i>						+	+
<i>S. enterica</i>						+	+
<i>S. marcescens</i>						+	+
<i>P. aeruginosa</i>						+	+
<i>A. baumannii</i>						+	+
<i>A. calcoaceticus</i>						+	+

Gram-Positive Bacteria	AMR Marker						
	SCCmec	mecA	mecC	vanA	vanB	CTX-M	KPC
Methicillin-resistant <i>Staphylococci</i>							
<i>S. aureus</i>	+	+	+				
<i>S. epidermidis</i>		+	+				
<i>S. lugdunensis</i>		+	+				
Vancomycin-resistant <i>Enterococci</i>							
<i>E. faecalis</i>				+	+		
<i>E. faecium</i>				+	+		

Plasma mNGS: Diagnostic Applications


Diagnostic Applications of Plasma mNGS




**Complicated
Pneumonia**




**Invasive Fungal
Infections**




FUO




Endocarditis




**Febrile
Neutropenia**



Osteomyelitis



Focal Abscess



Zoonoses

Diagnostic applications of plasma mNGS




Endocarditis

Eichenberger EM, et al.
(Duke University)

Fowler VG, et al.
(2023 Duke-ISCVID
Criteria)

Flurin L, et al.
(Mayo Clinic)




Pneumonia

Bergin SP, et al.
(PICKUP-
Multicentered)

Dworsky ZD, et al.
(Rady's Children)

Foong KS, et al.
(OSF St. Francis)




Invasive Fungal Infections

Armstrong AE, et al.
(Lurie Children's)

Vissichelli NC, et al.
(Virginia
Commonwealth
University)

Gracia, M., et al.
(Rady's Children)


Hoenigl M, et al



Febrile Neutropenia

Benamu E, et al.
(Stanford)

Yu et al.
(Advent Health Cancer
Institute)



Fever of Unknown Origin

Francisco DMA, et al.
(Baylor St. Luke)

Daher M, et al.
(Baylor St. Luke)


Shishido AA, et al.
(University of
Maryland Medical
Center)

Note that the provided information and clinical evidence is not exhaustive.

Considerations for Time of Use of the Plasma mNGS

Optimizing Plasma mNGS Utilization in Diagnostic Applications

CLINICAL TIMELINE



When to Consider Plasma mNGS

- Initial rapid testing is nondiagnostic
- Persistent or worsening infectious symptoms
- Contraindication to invasive diagnostic procedures
- Imaging evidence of infection
- High-risk for opportunistic or fastidious pathogens
- Critical illness with suspected infection


Results of Optimal Plasma mNGS Timing¹⁻⁵

- **Highest Pathogen DNA Concentration**
 - Collect samples early in infection or during critical illness for highest DNA concentration
- **Avoid Low Diagnostic Yield**
 - Prolonged, active antimicrobial therapy and/or source control can reduce diagnostic yield
- **Achieving Positive Clinical Impact**
 - Early testing has shown positive outcomes in retrospective analyses

1. Gaston et al. *Clin Lab Med.* 2024 Mar;44(1):63-73
2. Solanky et al. *Front. Trop. Dis.* 3:842100
3. Eichenberger et al. *Clin Infect Dis.* 2022 Jun 10;74(11):2020-2027
4. Shishido et al. *BMC Infect Dis.* 2022 Apr 13;22(1):372
5. Francisco et al. *Antimicrob Steward Healthc Epidemiol.* 2023 Feb 17;3(1):e31

Summary of Guidelines Including Plasma mNGS

Metagenomic Sequencing has Increasingly Been Incorporated into Major Infectious Disease Guidelines and Definitions




References:

1. Azar MM et al. *Am J Transplant*. 2022 Dec;22(12):3150-3169 2. Fowler VG et al. *Clin Infect Dis*. 2023 Aug 22;77(4):518-526 3. DiSimone DC et al. *J Am Heart Assoc*. 2025;14:e040218 4. Douglas AP et al. *Transplant Cell Ther*. 2025;31(4):194-223 5. Panel on Guidelines for the Prevention and Treatment of Opportunistic Infections in Adults and Adolescents With HIV. National Institutes of Health, HIV Medicine Association, and Infectious Diseases Society of America. Accessed June 6th, 2025 6. Prevention and Treatment of Cancer-Related Infections. Version 1.2025, 06/20/2025 © 2025 National Comprehensive Cancer Network® (NCCN®)


BALF Metagenomic Sequencing Product Overview

BALF Metagenomic Sequencing (mNGS) of microbial cell-free DNA (mcfDNA)



- Bronchoalveolar lavage fluid-based mcfDNA metagenomic sequencing test
- Designed for broad detection and classification of over 500 pathogens associated with lung infections

Patients who may be appropriate for **BALF mNGS**




- Suspected lung infection
- Immunocompromised
- Undergoing bronchoalveolar lavage (BAL)

mNGS's broad approach is optimized for adult & pediatric **immunocompromised** patients, are more likely to harbor difficult-to-diagnose pathogens.

Test Process Overview

Reports are typically available one day after specimen receipt



BALF mNGS can rapidly identify the potential etiology of pneumonia and improve diagnostic yield over SOC testing



Harnessing
mcfDNA from a
low volume of BAL
fluid sample
(1 mL)

Convenient



>50% increased
relative diagnostic
yield compared to
all SOC-BAL tests
combined¹

Comprehensive



Rapid identification of
pathogens with
results typically
available within 1 day
of sample receipt

Timely

1. Unpublished internal data from analysis performed on patients enrolled in PICKUP study; n=118; BALF mNGS identified adjudicated cause of pneumonia in 40 patients versus 29 with SOC alone

Questions?