



Department of
Health

Clinical Microbiology

A Refresher

Objectives


- Explain the infectious diseases diagnostic testing process, including appropriate specimen collection.
- Describe the bacterial/fungal culture and antimicrobial susceptibility process.
- Define beta-lactamases and discuss how they contribute to bacterial resistance.

Outline

- Basics of microbiology laboratory testing
- Bacterial/fungal cultures
- Antimicrobial susceptibility testing
- Antimicrobial resistance



Basics of Microbiology Laboratory Testing



Specimen Collection

- Most important part of the testing processing
- Good test results begin with a good specimen
- Includes:
 - Collection of specimen
 - Storage of specimen
 - Transport of specimen

Specimen Collection: Appropriate Specimen

- The specimen collected from a patient depends on clinical symptoms
 - Blood
 - Cerebrospinal fluid (CSF)
 - Respiratory
 - Throat
 - Sputum
 - Bronchial washings
 - Eye
 - Ear
 - Wound / Tissues
 - Urine
 - Genital specimens
 - Stool

Sterile vs. non-sterile site

Specimen Collection: Important Factors

- Timing of collection
- Amount of material
- Contamination
 - Normal microbiota or other contamination
- Labeling
- Storage
- Transportation

Case in Point: *C. difficile*

- Collect specimens only on symptomatic patients
 - Colonization vs. Infection
- Collect stool
 - Unformed, liquid stool
 - Patient should not be on laxatives
 - Collect in approved sterile container based on lab guidance
 - Transport to lab within designated time and temperature

Fig 1. **Bristol Stool Chart**

Type 1		Separate hard lumps, like nuts (hard to pass)
Type 2		Sausage-shaped but lumpy
Type 3		Like a sausage but with cracks on its surface
Type 4		Like a sausage or snake, smooth and soft
Type 5		Soft blobs with clear cut edges (passed easily)
Type 6		Fluffy pieces with ragged edges, a mushy stool
Type 7		Watery, no solid pieces. Entirely liquid

Infectious Diseases Diagnostic Tests

- Microscopy
- Culture
- Molecular methods
- Immunological analysis

Microscopy

- Can be performed directly from patient specimen or from a culture
- Mainly used for bacteria, fungi and parasites
 - Rarely used for viruses
- Various stains may be used to visualize the microorganism
 - Example: Gram stain

Culture

- Mainly used for bacteria, fungi
- Rarely used for viruses
 - Viruses are obligate intracellular pathogens
 - Need host cells to grow; must grow on live cells in the laboratory
- Some parasite can be cultured but it's rarely performed

Molecular Methods

- Can be performed directly from patient specimen or from a culture
- Used for bacteria, fungi, viruses and parasites
- Includes:
 - DNA probes
 - NAAT (nucleic acid amplification test)
 - PCR (polymerase chain reaction)
 - 16s
 - Sequencing

Immunological Analysis

- Evaluation of
 - Antigen - protein on microorganism
 - Antibody - immune response to microorganism
 - IgM antibody – produced during acute phase of disease
 - IgG antibody – produced later on during the disease process; can persist for life
- Examples test methods
 - ELISA, EIA, lateral flow

Microscopy

- Can be performed directly from patient specimen or from a culture
- Mainly used for bacteria, fungi and parasites
 - Rarely used for viruses
- Various stains may be used to visualize the microorganism
 - Example: Gram stain

Culture

- Mainly used for bacteria, fungi
- Rarely used for viruses
 - Viruses are obligate intracellular pathogens
 - Need host cells to grow; must grow on live cells in the laboratory
- Some parasite can be cultured but it's rarely performed

Molecular Methods

- Can be performed directly from patient specimen or from a culture
- Used for bacteria, fungi, viruses and parasites
- Includes:
 - DNA probes
 - PCR (polymerase chain reaction)
 - NAAT (nucleic acid amplification test)
 - 16s
 - Sequencing

Immunological Analysis

- Evaluation of
 - Antigen - protein on microorganism
 - Antibody - immune response to microorganism
 - IgM – acute response
 - IgG – chronic or resolved infection
- Used for bacteria, fungi, viruses and parasites




Bacterial/Fungal Culture

Culture

- The process of growing microorganisms by taking bacteria from the infection site and growing them in the artificial environment of the laboratory
- Must recreate body environment in the lab including
 - Temperature
 - Atmosphere
 - Nutrients

Culture Media

- Provide the nutrients the organism needs to grow outside the body



Harriott, 2012

Culture Process


- Process patient specimen
- Choose correct media plates for specimen type and the organism that may be causing the infection
- Inoculate media with patient specimen
- Incubate at correct temperature, atmosphere and for the correct amount of time

Culture Incubation Time


- Appropriate time
 - Varies for specimen and organism
 - 1 -2 days for most bacteria
 - 3 weeks for fungi
 - 6 weeks for mycobacteria
- After predetermined incubation time, remove plates and identify any organisms growing

Identification of Growth

- Evaluate plates after designated incubation period for growth of pathogens
- Must distinguish pathogen from microbiota




Harriott, 2013



Harriott, 2013

Normal Flora



<https://microbenotes.com/normal-human-microbiota/>

Identification of Growth

- If a potential pathogen is suspected, then further testing is warranted
- Type of tests
 - Phenotypic
 - Evaluation of the expression of genetic characteristics
 - Genotypic
 - Use of molecular techniques to evaluate RNA or DNA from colony growing on plate


Phenotypic Identification

- Evaluation of colonies
- Biochemical testing
 - Mostly automated
- MALDI-TOF
 - Mass spectroscopy

Staphylococcus aureus



Harriott, 2012





Harriott, 2012




Harriott, 2004

MALDI-TOF






Blood Cultures



Culture Reports

- Sterile sites
 - Growth or no growth
- Normal flora
- Quantity of organisms
 - Rare, Few, Moderate, Many
 - 1+, 2+, 3+, 4+
 - Urines: colony forming units per mL of urine (CFU/mL)

Growth Quantitation





Antimicrobial Susceptibility Testing and Stewardship

Antimicrobial Susceptibility Testing

- Testing to determine if an antimicrobial agent can be used to effectively treat an infection
- Phenotypic
 - Evaluating expression of a gene
 - Main method for bacteria and fungi
- Genotypic
 - Evaluating the presence of a gene
 - Can be performed for bacteria, fungi, and viruses
- Susceptibility testing for parasites rarely performed

Phenotypic Susceptibility Testing

- Purpose:
 - Predict whether the bacteria is capable of expressing resistance to the antimicrobial agents that are potential choices as therapeutic agents for managing the infection
- The tests are performed under standardized conditions
 - Ensures reproducibility
- Guide not a standard!!
 - Results should be combined with clinical information and experience when selecting the most appropriate antibiotic for a particular patient.

Choosing Antimicrobials to Test

- Number of factors play a role in this
 - Gram-negative vs. Gram-positive
 - Intrinsic resistance
 - Availability of antimicrobials
 - Site of infection

Main Methods

- Disk diffusion (Kirby-Bauer)
 - Broth Dilution
 - E-test
-
- These methods are standardized to ensure reproducibility and accuracy


Disk Diffusion (Kirby Bauer)

- Concentrations of antimicrobial agents are incorporated into filter paper disk
- Culture plate is incubated with bacteria and discs overnight
- Measure zone of inhibition and determine if it is susceptible, intermediate, or resistant



Broth Dilution Method

- Different concentrations of antibiotics are placed in micro dilution tray
- A suspension of the bacterial organism is also added to each wells
- Incubate overnight
- Look for growth or lack of growth in well (turbidity indicates growth)
- Determine the MIC






Antimicrobial Resistance

Antibiotic Resistance Mechanisms

- Altered antibiotic transport
 - Efflux pumps
- *Enzyme inactivation*
 - *Activity of the antibiotic is hindered due to hydrolysis*
- Target by-pass
 - Prevent the drug from getting to the target site
- Target modification
 - Bacteria modify the antibiotic target rendering it inactive
- These mechanisms can be
 - inherited (passed on generation to generation)
 - acquired




Bacterial Genome

- Chromosomal DNA
- Extrachromosomal DNA (Plasmid)
 - Physically separated from chromosome
 - Contains few genes, not essential for the organism to function
 - Replicates independently of the genome
 - Type of mobile genetic element
 - Can be horizontally transferred


Production of Beta-Lactamase Enzymes

- Bacteria destroy the beta-lactam ring
 - The antibiotic is no longer active
- Examples:
 - Extended spectrum beta-lactamase (ESBL)
 - AmpC beta-lactamase
 - Carbapenemase

β -Lactam Antibiotics



Penicillins



Cephalosporins

Rocha, Antônio & Barsottini, Mario & Rocha, Renan & Laurindo, Maria & Leandro, Francisco & Moraes, Francisco & Rocha, Soraya. (2019). *Pseudomonas Aeruginosa: Virulence Factors and Antibiotic Resistance Genes*. Brazilian Archives of Biology and Technology. 62. 19180503. 10.1590/1678-4324-2019180503.

Carbapenemase Genes

The “Big-5”

- *Klebsiella pneumoniae* carbapenemase (KPC)
- Imipenemase Metallo- β -lactamase (IMP)
- New Delhi Metallo- β -lactamase (NDM)
- Oxacillinase-48-type carbapenemase (OXA-48)
- Verona integron-encoded Metallo- β -lactamase (VIM)

Summary

- Identification of microorganisms from patient specimen begins in specimen collection
 - Most important step
- Types of tests include culture, molecular and immunological
- Antimicrobial susceptibility testing is performed to help determine which antimicrobial can be used for treatment
 - Testing can be phenotypic or genotypic
- Beta-lactamases are enzymes produced by some bacteria that inactivate antibiotics resulting in resistance



Questions?
Melphine.Harriott@tn.gov